About us

ZOE BIO

Elevating Lives through **Exceptional Biotechnology**

At ZOE BIO, our driving force is our unwavering commitment to harnessing cutting-edge biotechnology to create a more meaningful and valuable life for our consumers.

We push the boundaries of bioengineering to empower individuals to lead more fulfilling lives, believing that true beauty goes beyond skin-deep. Our advanced biomedical technology not only restores beauty but also promotes holistic well-being.

With ZOE BIO, you can trust that we are devoted to helping you achieve a more beautiful and meaningful life.

Experience the transformative path to redefining beauty and shaping a brighter future.



formule

Prescription beyond beauty



Treatment application area

fomule F2 • fomule F3





ZOE BIO Inc.

1F, Moohan Building, 54, Bongeunsa-ro 20-gil, Gangnam-gu,

Z = BIO

www.zoebio.co.kr

HELP Technology

HA Dermal Filler

formule F1 F2 F3

High molecular weight HA **E**nhanced visco-elasticity Low crosslinking agent Pure gel

The HELP Technology is to maximize the safety and effectiveness of HA filler and to minimize the change form HA origin.



Superficial dermis	Deep dermis	
Soft	Medium	Hard

High Cross-linking Efficacy Technology

Use of small amount of cross-linking agent & High viscoelasticity → Increased safety & effectiveness

Pure Gel Technology

Removal of residual BDDE and impurities & swelling → Low cross linking agent with high safety

Homogeneous Gel Technology

Easy to control with stable extrusion force

formule P/HP/PCT

formule P

Pure PN only

- Stimulation of fibroblast and metabolic activity
- Tissue regeneration efficacy
- Stimulation of the cell production
- Stimulation of the secretion of various growth factors

formule HP

PN + HA collaboration

Increase of HA functions in moisturizing & adding volume in the skin



formule PCT

High regenerative support

skin booster with safety and easy to use



- Benefits Improves skin elasticity
 - Improves moisturizing
 - Reduces fine lines and wrinkles
 - Improves skin barrier function
 - Helps skin normalizing
 - Improves skin tone & texture


formule **HiLo**

High&Low Molecule Skin Booster

formule Hilo

80mg HA

- 40 mg of High MW HA
- 40 mg of Low MW HA



- Highly purified HA (safety)
- Visco-supplement
- Low viscosity / Optimal tissue diffusion
- Long duration as comparable with low cross-linked gel
- Bio-supplement
- Controlled inflammatory response
- Activation of cell-remodeling

Benefits

- Activates multi-level dynamic remodeling and improves skin elasticity of face, neck and body.
- Improves the extracellular environment (in vitro studies)

Products Specification

HA Dermal Filler				
Classification	formule F1	formule F2	formule F3	
Ingredient	HA 24 mg/ml Lidocaine 0.3 %			
Application	Superficial Dermis	Mid/Deep Dermis	Deep Dermis	
Volume		1.1 ml		
Needle	30G	27G	25G	
Storage Condition		2 ~ 25 °C		
Expiration Date		24 months		
Package Unit	1 Syrir	nge + 2 Needles / 1 I	Вох	

Polynucleotide(PN) Skin Booster				
Classification	formule P	formule HP	formule PCT	
Ingredient	PN 20 mg/ml	HA 10 mg/ml PN 10 mg/ml Lidocaine 0.3 %	Hyaluronic acid + Polynucleotid	
Application	Superficial Dermis			
Volume	2.2 ml	2.2 ml	2.5 ml	
Needle	Not included (33G Recommended)		Not included (31G Recommende	
Storage Condition	2 ~ 25 °C			
Expiration Date	24 months			
Package Unit	2 Syringes /1 Box	1 Syringe / 1 Box	3 Syringes / 1 Box	

High & Low Molecule HA Skin Booster

Classification	tormule HiLo	
Ingredient	HA 80 mg/2.5 ml/Syr. [HW-HA 40mg + LW – HA 40mg]	
Application	Superficial Dermis	
Volume	2.5 ml	
Needle	Not included (29G Recommended)	
Storage Condition	2 ~ 25 °C	
Expiration Date	24 months	
Package Unit	1 Syringe / 1 Box	